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1. Ritz J, Pesando JM, Notis-McConarty J, Lazarus H, Schlossman SF: 1980. A monoclonal antibody to human acute lymphoblastic leukemia antigen. Nature 283 : 583.
2. Ritz J, Nadler LM, Bhan AJ, Notis-McConarty J, Pesando JM, Schlossman SF: 1981 . Expression of common acute lymphoblastic leukemia antigen (CALLA) by lymphomas of B-cell and T-cell lineage. Blood 48: 648.
3. Clavell LA, Lipton JM, Blast RC. Kudish M, Pesando JM, Schlossman SF, Ritz J: 1981. Absence of common ALL antigen on normal bipotent myeloid, erythroid and granulocytic progenitors. Blood 58 : 333.
4. Sallan SE, Ritz J, Pesando JM, Gelber R. O'Brien C, Hitchcock S, Coral FS, Schlossman SF: 1980. Cell surface antigens: Prognostic implication in childhood acute lymphoblastic leukemia. Blood 55 : 395.
5. Greaves MF: 1981. Monoclonal antibodies as probes for leukemic heterogeneity and hemopoetic differentiation. In: Leukemia Markers, ed. W. Knapp, Academic Press, p. 19.
6. Greaves MF: 1981. Biological basis of differential prognosis in immunologically defined sub-classes of acute lymphoblastic leukemia. In: Leukemia markers, ed. W. Knapp, Academic Press, p. 463.
7. Nadler LM, Ritz J, Todd RFIII. Griffin J, Reinherz EL, Schlossman SF: 1981. Diagnosis and treatment of human leukemia and lymphomas utilizing monoclonal antibodies. Progress in Hematology, ed. E. Brown, Grune and Stratton, Inc., NY., p. 187.
8. Cossman J, Neckers LM, Leonard WJ, Greene WC: 1983. Polymorphonuclear neutrophils express the common acute lymphoblastic leukemia antigen. J Exp Med 157 : 1064.
9. Hokland P, Ritz J, Schlossman SF, Nadler LM: 1985. Orderly expression of B cell antigens during the in vitro differentiation of nonmalignant human pre-B cells. J Immunol 135 : 1746.
10. Boucheix C, Perrot JY, Mirshahi M, Giannoni M, Billard F, Bernadou A, Rosenfeld C: 1985. A new set of monoclonal antibodies against acute lymphoblastic leukemia. Leuk res 9(5) : 597-604.
11. Boucheix C. Perrot JY, Mirshahi M, Bernadou A. Rosenfeld C: 1983. A rapid method for detection of membrane amtigens by Fluorescent microscopy and its application to screening hybridoma antibodies. J Immunol Methods 57 : 145-150.
12. Lebien WT, Mac. Cormach RT: 1989. The common acute lymphoblastic leukemia antigen (CD10). Emancipation from a functional Enigma. Blood 73(3) : 625-635.
13. Penchansky L, Pirrotta V, Kaplan SS: 1993. Flow cytometric study of the expression of neutral endopeptidase (CD10/CALLA) on the surface of newborn granulocytes. Mod Pathol 6:414-418.
14. Connelly JC. Chambless R, Holiday D, Chittenden K, Johnson AR : 1993. Up-regulation of neutral endopeptidase (CALLA) in human neutrophils by granulocyte-macrophage colony-

stimulating factor. J Leukoc Biol 53:685-690.
15. Taki T, Ida K, Bessho F, Hanada R, Kikuchi A, Yamamoto K, Sako M, Tsuchida M, Seto M, Ueda R, Hayashi Y: 1996. Frequency and clinical significance of the MLL gene rearrangemants in infant acute leukemia. Leukemia 10:1303-1307.
16. Heerema NA, Maben KD, Bernstein J. Breitfeld PP, Neiman RS, Vance GH: 1996. Dicentric (9;20)(p11;q11) identified by fluorescence in situ hybridization in four pediatric acute lymphoblastic leukemia patients. Cancer Genet Cytogenet 92:1 11-115.
17. Ganju RK. Shpektor RG, Brenner DG, Shipp MA:1996. CD10/neutral endopeptidase 24. 11 is phosphorylated by casein kinase ll and coassociates other phosphoproteins including the lyn src-related kinase. Blood 88:4159-4165.

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