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1. Centers for Disease Control and Prevention. *Acute Traumatic Spinal Cord Injury Surveillance -- United States, 1987.* http://www.cdc.gov/mmwr/preview/mmwrhtml/00000022.htm
2. National Spinal Cord Injury Statistical Center. *THE 2004 ANNUAL STATISTICAL REPORT for the MODEL SPINAL CORD INJURY CARE SYSTEMS.* http://www.spinalcord.uab.edu
3. Bravo, G., Guizar-Sahagun, G., Antonio, I., David. C., Carlos M. (2004). *Cardiovascular alternations After Spinal Cord Injury: An overview. Curr. Med. Chem-Cardiovascular & Hematological agents.* 2: 133-148.
4. John F. Bilello, MD. (2003) *Cervical spinal Cord Injury and the Need for Cardiovascular Intervention.* American Medical Association. 10: 1127-1129.
5. Irene S. Gilgoff, MD. (1991) *Cardiac Pacemaker in High Spinal Cord Injury.* Arch Phys Med Rehabll Vol 72. 72: 601-3.
6. Devivo, M. J., Go, B. K., & Jackson, A. B. (2002). *Overview of the national spinal cord injury statistical center database.* Journal of Spinal cord Medicine, 25(4): 335-338.
7. Hoitt, J. D., Banzett, R. B., Lohmeir, H. L., Hixon, T. J., & Brown, R. (2003). *Clinical ventilator adjustments that improve speech.* Chest, 124(4): 1512-1521.
8. Krause, J. S., Devivo, M. J., & Jackson, A. B. (2004). *Health status, community integration, and economic risk factors for mortality after spinal cord injury.* Archives of Physical Medicine and Rehabilitation, 85(11): 1764-1773.
9. Carter, R. E., Donovan, W. H., Halstead, L., & Wilkerson, M. A. (1987). *Comparative study ofelectophrenic nerve stimulation and mechanical ventilatory support in traumatic spinal cordinjury.* Paraplegia, 25: 86-91.
10. Wolf, C., Meiners, T., & Eisenhuth, J. (1999, 2000). *Personal situation and situation at home of respiratory-dependent patients with tetraplegia: electrophrenic versus mechanical ventilation.* Paper presented at the Functional Electrical Stimulation in Paralysed Respiratory Muscles Proceedings of the Workshop, Hamburg Germany
11. Higgins, P. A. (1998). *Patient perception of fatigue while undergoing long-term mechanical ventilation: Incidence and associated factors.* Heart & Lung, 27(3): 177-183.
12. Rossi Ferrario, S., Zotti, A. M., Zaccaria, S., & Donner, C. F. (2001). *Caregiver stain associated with tracheostomy in chronic respiratory failure.* Chest. 119(5): 1498-1502.
13. van Kesteren, R. G., Velthuis, B., & van Leydan, L. W. (2001). *Psychosocial problems arising from home ventilation.* American Journal of Medical Rehabilitation, 80(6): 439-496.
14. Anderson, Kim (2004). *Targeting Recovery: Priorities of the spinal cord-Injured Populaton.* Journal of Neurotrauma 21(10): 1371-1383.
15. William W. L. Glenn, M.D. (1972) *Total Ventilatory support in a quadriplegic patient with radiofrequency electrophrenic respiration.* 286(10): 513-516.
16. Dimarco AF, Onders RP, Kowalski KE, Miller ME, Ferek S, Mortimer JT. (2002) *Phrenic nerve pacing in a tetraplegic Patient via intramuscular diaphragm electrodes.* American Journal of Respiratory and Critical Care Medicine. 166: 1604-1606.
17. Glenn, W. W., J. F. Hogan, et al. (1980). *Ventilatory support of the quadriplegic patient with respiratory paralysis by diaphragm pacing.* Surg-Clin-North-Am 60(5): 1055-78. issn: 0039-6109.
18. Anthony F. DiMarco. (1999). *Diaphragm Pacing in Patients with Spinal Cord Injury.* Top Spinal Cord Inj Rehabil. 5(1): 6-20.
19. John Moxham., John M. Shneerson. (1993) *Diaphragmatic Pacing.* Am Rev Respir Dis. 148: 533-536.
20. Glenn, W. W., J. F. Hogan, et al. (1984). *Ventilatory support by pacing of the conditioned diaphragm in quadriplegia.* N-Engl-J-Med. 310(18): 1150-5 issn: 0028-4793.
21. AH Moss. (1996) *Patients with amyotrophic lateral sclerosis receiving long-term medical ventilation. Advance care planning and outcomes.* CHEST. The Cardiopulmonary and Critical Care Journal. 110: 249-255.
22. Creasey, G., J. Elefteriades, et al. (1996). *Electrical stimulation to restore respiration.* J-Rehabil-Res-Dev 33(2): 123-32 issn: 0748-7711.
23. Nochomovitz, M. L., D. K. Peterson, et al. (1988). *Electrical activation of the diaphragm.* Clin-Chest-Med 9(2): 349-58 issn: 0272-5231.
24. Glenn, W. W., R. T. Brouillette, et al. (1988). *Fundamental considerations in pacing of the diaphragm for chronic ventilatory insufficiency: a multi-center study.* Pacing-Clin-Electrophysiol 11(11 Pt 2): 2121-7 issn: 0147-8389.
25. Synapse Biomedical Inc. "Memo" (2007)
26. Aiyar, H., T. A. Stellato, et al. (1999). *Laparoscopic Implant Instrument for the Placement of Intramuscular Electrodes in the Diaphragm.* IEEE Trans Rehab eng 7(3): 360-71.
27. David K. Peterson. (1994) *Long-Term Intramuscular Electrical Activation of the Phrenic Nerve: Safety and Reliability.* IEEE Transactions of Biomedical engineering. 41(12): 1115-26.
28. Schmit, B. D. and J. T. Mortimer (1993). *Electrical Activation of the Diaphragm using Epimysial Electrodes.* 15th Int. Conf. IEEE-EMBS. 1251-2.
29. Schmit, B. D., T. A. Stellato, et al. (1998). *Laparoscopic Placement of Electrodes for Diaphragm Pacing Using Stimulation to Locate the Phrenic Nerve Motor Points.* IEEE Trans Rehab Eng 6(4): 382-390.
30. Onders RP., Aiyar H., Mortimer JT. (2004) *Characterization of the human diaphragm muscle with respect*

to the phrenic nerve moter points for diaphragmatic pacing. Am. Surg. **70**(3):241-7.

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1. David K. Peterson, et al., (1994) *Long-Term Intramuscular Electrical Activation of the Phrenic Nerve: Efficacy as a Ventilatory Prosthesis.* IEEE Transactions on Biomedical Engineering.**41**(12): 1127-35.
2. David K. Peterson, et al., (1994) *Long-Term Intramuscular Electrical Activation of the Phrenic Nerve: Safety and Reliability.* IEEE Transactions on Biomedical Engineering. **41**(12): 1115-26.
3. E. IN VIVO EVALUATION/E.2/3 Movement and Recrutiment Stability[redacted]
4. E. IN VIVO EVALUATION/E.2/3 Movement and Recrutiment Stability[redacted]

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


[redacted]

[redacted]

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1. Meyer E., Sohr D., Gastmeier P. et al., (2009) *New identification of outliers and ventilator-associated pneumonia rates from 2005 to 2007 within the German Nosocomial Infection Surveillance System.* J. Hosp. Infect. **73**(3): 246-52.
2. Margaret AD., Lindsey M., Weiner MPH., et al., (2013) *National Healthcare Safety Network (NHSN) report, data summary for 2012, Device-associated module.* Am. J. Infect. Control. **41**: 1148-66.
3. The Acute Respiratory Distress Syndrome Network, (2000) *Ventilation with lower tidal volume as compared with traditional tidal volumes for acute lung injury and the acute respiratory distress syndrome.* New Eng. J. Med. **342**(18): 1301-1308.

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国　名	販売名	使用目的	承認年月日
米国(FDA)	NeuRx DPS RA/4 Diaphragm Pacing Stimulation System	(SCI) This device is indicated for use in patients with stable, high spinal cord injuries with stimulatable diaphragms, but lack control of their diaphragms.	2008/6/17
	NeuRx DPS Diaphragm Pacing System	(ALS) This device is indicated for use in amyotrophic lateral sclerosis (ALS) patients with a stimulatable diaphragm (both right and left portions) as'demonstrated by voluntary contraction or phrenic nerve conduction studies, and who are experiencing chronic hypoventilation (CH), but not progressed to an FVC less than 45% predicted.	2011/9/28
欧州(CE) ・ドイツ ・フランス ・イギリス ・オランダ ・スペイン ・ベルギー ・イタリア ・スイス ・ノルウェー ・アイスランド	The NeuRx Diaphragm Pacing System™ (NeuRx DPS™)	(横隔膜機能不全) The NeuRx Diaphragm Pacing System™ (NeuRx DPS™) provides ventilatory support in patients with diaphragm dysfunction of neuromuscular origin. Diaphragm dysfunction can result in abnormal or absent respiration in patient populations of high-level spinal cord injury, central hypoventilation, amyotrophic lateral sclerosis, and other injuries or diseases affecting the neuromuscular respiratory pathways.	2007/9/5
サウジアラビア	NeuRx™ diaphragm pacing system. (Intramuscular percutaneous motor point diaphragm stimulation system)	(横隔膜機能不全) The NeuRx Diaphragm Pacing System™ (NeuRx DPS™) provides ventilatory support in patients with diaphragm dysfunction of neuromuscular origin. Diaphragm dysfunction can result in abnormal or absent respiration in patient populations of high-level spinal cord injury, central hypoventilation, amyotrophic lateral sclerosis, and other injuries or diseases affecting the neuromuscular respiratory pathways.	2010
カナダ	NEURX Diaphragm Pacing System – DPS System Kit	(ＳＣＩ) The NeuRx Diaphragm Pacing System™ (NeuRx DPS™) provides ventilatory support in patients with diaphragm dysfunction of neuromuscular origin. Diaphragm dysfunction can result in abnormal or absent respiration in patient populations of high-level spinal cord injury, central hypoventilation, and other injuries or diseases affecting the neuromuscular respiratory pathways.	2013/3/20
オーストラリア	NIntramuscular diaphragm/phrenic nerve electrical stimulation system	(横隔膜機能不全) The NeuRx Diaphragm Pacing System is intended to provide ventilatory support in patients with diaphragm dysfunction of nueromuscular origin and is implanted in an outpatient laparoscopic procedure. The NeuRX Diaphragm Pacing System helps patients breathe by stimulating their diaphragm muscles, through repetitive electrical	2009/1/20

		stimulation to the implanted electrodes to cause the patient's diaphragm to contract and cause the patient to draw breath in a manner similar to natural breathing.	
イスラエル	NeuRx™ diaphragm pacing system. (Intramuscular percutaneous motor point diaphragm stimulation system)	(横隔膜機能不全) The NeuRx Diaphragm Pacing System™ (NeuRx DPS™) provides ventilatory support in patients with diaphragm dysfunction of neuromuscular origin. Diaphragm dysfunction can result in abnormal or absent respiration in patient populations of high-level spinal cord injury, central hypoventilation, amyotrophic lateral sclerosis, and other injuries or diseases affecting the neuromuscular respiratory pathways.	2008/ 11/3
コロンビア	NeuRx™ diaphragm pacing system. (Intramuscular percutaneous motor point diaphragm stimulation system)	(横隔膜機能不全) The NeuRx Diaphragm Pacing System™ (NeuRx DPS™) provides ventilatory support in patients with diaphragm dysfunction of neuromuscular origin. Diaphragm dysfunction can result in abnormal or absent respiration in patient populations of high-level spinal cord injury, central hypoventilation, amyotrophic lateral sclerosis, and other injuries or diseases affecting the neuromuscular respiratory pathways.	2013/ 7/11
トルコ	NeuRx™ diaphragm pacing system. (Intramuscular percutaneous motor point diaphragm stimulation system)	(横隔膜機能不全) The NeuRx Diaphragm Pacing System™ (NeuRx DPS™) provides ventilatory support in patients with diaphragm dysfunction of neuromuscular origin. Diaphragm dysfunction can result in abnormal or absent respiration in patient populations of high-level spinal cord injury, central hypoventilation, amyotrophic lateral sclerosis, and other injuries or diseases affecting the neuromuscular respiratory pathways.	2012
ブラジル	NeuRx™ diaphragm pacing system. (Intramuscular percutaneous motor point diaphragm stimulation system)	(ALS) Treatment of ALS	2010/ 4/6
ヨルダン	NeuRx Pacing System	(横隔膜機能不全) The NeuRx Diaphragm Pacing System™ (NeuRx DPS™) provides ventilatory support in patients with diaphragm dysfunction of neuromuscular origin. Diaphragm dysfunction can result in abnormal or absent respiration in patient populations of high-level spinal cord injury, central hypoventilation, amyotrophic lateral sclerosis, and other injuries or diseases affecting the neuromuscular respiratory pathways.	2011/ 10/27

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