Pharmaceuticals and Medical Devices Agency Office of Review Management December 1, 2025

Posting Date	Monograph	Test	Brand Name
December 1, 2025	Solifenacin Succinate	Purity (1) Related substances 1, Purity (2) Related substances 2, Purity (3) Related substances 3, Assay	Develosil ODS-UG-5
		Purity (5) Enantiomer and diastereomer	CHIRALPAK AD-H
December 1, 2025	Solifenacin Succinate Orally Disintegrating Tablets	Purity Related substances, Uniformity of Dosage Unit, Dissolution, Assay	Develosil ODS-UG-5
December 1, 2025	Solifenacin Succinate Tablets	Purity Related substances, Uniformity of Dosage Unit, Dissolution, Assay	Develosil ODS-UG-5
September 1, 2025	Gliclazide Tablets	Uniformity of dosage unit, Assay	NUCLEOSIL 120-5 C18
June 2, 2025	Dilazep Hydrochloride Tablets	Identification, Uniformity of dosage units, Assay	Develosil ODS-7
December 2, 2024	Tandospirone Citrate	Purity	Inertsil ODS-3(Model No.5020-01771)
December 2, 2024	Tandospirone Citrate Tablets	Dissolution, Assay	L-columm ODS (Model No.622070)
September 30, 2024	Aspirin Enteric-coated Tablets	Purity, Uniformity of dosage units, Assay	LiChrospher 100 RP-18 KANTO CHEMICAL (Ser.No.50823-1B) Merck KGaA (CAT No.1.50823.0001)
September 2, 2024	Aripiprazole Tablets	Identification, Uniformity of dosage units, Assay	YMC-Pack ODS-A or YMC-Pack ODS- AM
		Purity Related substances Dissolution	Develosil ODS-HG-5 YMC-Pack ODS-A

Pharmaceuticals and Medical Devices Agency Office of Review Management December 1, 2025

Posting Date	Monograph	Test	Brand Name
September 2, 2024	Blonanserin	Purity Related substances, Assay	L-column ODS (Model No.622070)
September 2, 2024	Blonanserin Powder	Uniformity of dosage units, Dissolution, Assay	L-column ODS(Model No.622070)
September 2, 2024	Blonanserin Tablets	Uniformity of dosage units, Dissolution, Assay	L-column ODS(Model No.622070)
September 2, 2024	Celecoxib Tablets	Uniformity of dosage unit, Assay	SUPELCOSIL LC-DP,5µm, 4.6 × 250mm
September 2, 2024	Dapagliflozin Propylene Glycolate	Identification (2), Purity Related substances, Assay	Waters SunFire C18
September 2, 2024	Hydrate	Propylene Glycol	Restek / Carbowax, Stabilwax Crossbonded
September 2, 2024	Dapagliflozin Propylene Glycolate Tablets	Identification (2), Purity Related substances, Uniformity of dosage units, Assay	YMC-Pack Pro C18
		Dissolution	Phenomenex Luna C18(2)
		Purity (3) Related substances	Novapak
September 2, 2024	Esomeprazole Magnesium Hydrate	Enantiomer	Merck Chiral AGP
		Assay	Merck LiChrophere
		Identification	Merck Chiral AGP
September 2, 2024	Esomeprazole Magnesium Delayed-	Purity	Microspher C18
	Release Capsules	Uniformity of dosage units, Dissolution, Assay	Zorbax Eclipse XDB C18
September 2, 2024	Limaprost Alfadex Tablets	Purity Related substances, Uniformity of dosage units, Assay	YMC-Pack ODS-A or Develosil ODS-5

Pharmaceuticals and Medical Devices Agency Office of Review Management December 1, 2025

Posting Date	Monograph	Test	Brand Name
September 2, 2024	Perospirone Hydrochloride Hydrate	Purity Related substances, Assay	YMC-Pack ODS AM (Model No.AM12S05-1506WT)
September 2, 2024	Perospirone Hydrochloride Tablets	Uniformity of dosage units, Dissolution, Assay	YMC-Pack ODS AM (Model No.AM12S05-1506WT)
September 2, 2024	Pregabalin	Purity (1) Related substances (i) Polar related substances, Purity (1) Related substances (ii) Nonpolar related substances, Assay	YMC-Pack ODS-AQ 5μm 4.6 × 250 mm
		Purity (2) Enantiomer	Hypersil BDS C18 5µm 4.6 × 250 mm
September 2, 2024	Pregabalin Capsules	Purity Related substances, Uniformity of dosage units, Assay	μ-Bondapak C18 10μm 3.9 × 300 mm
		Dissolution	ZORBAX SB-CN 5µm 4.6 × 150 mm
September 2, 2024	Pregabalin Orally Disintegrating Tablets	Purity Related substances, Uniformity of dosage units, Assay	XBridge BEH Shield RP18 3.5µm 4.6 × 100 mm
		Dissolution	ZORBAX SB-CN 5µm 4.6 × 150 mm
June 3, 2024	L-Carbocisteine for Syrup	Uniformity of dosage units, Dissolution, Assay	Inertsil ODS-2
June 3, 2024	Levetiracetam	Purity (2) Related substances, Assay	ChromSpher Si Luna Silica(2)
		Purity (3) Enantiomer	CHIRALCEL OD
		Identification (2), Purity (1) Related	Symmetry C18
June 3, 2024	Levetiracetam tablets	substances, Assay	Symmetry Shield C18
		Purity (2) Enantiomer	CHIRALCEL OD

Pharmaceuticals and Medical Devices Agency Office of Review Management December 1, 2025

Posting Date	Monograph	Test	Brand Name
June 3, 2024	Midodrine Hydrochloride	Purity (2) Related substances, Assay	COSMOSIL 5C ₁₈ -AR-II
	Midodrine Hydrochloride Orally	Purity Related substances	Develosil TMS-UG-5
June 3, 2024	Disintegrating Tablets	Uniformity of dosage units, Dissolution, Assay	TSKgel ODS-80Ts
		Purity Related substances	Develosil TMS-UG-5
June 3, 2024	Midodrine Hydrochloride Tablets	Uniformity of dosage units, Dissolution, Assay	TSKgel ODS-80Ts
Danamban 1 0000	Compostet Medilete Tablete	Uniformity of dosage units	Finepak SIL C18S
December 1, 2023	Camostat Mesilate Tablets	Assay	Finepak SIL C18S
Danamilan 1 2002	Pranlukast Capsules	Uniformity of dosage units	Develosil C8-HG-5
December 1, 2023		Assay	Develosil C8-HG-5
December 1, 2023	Pranlukast for Syrup	Uniformity of dosage units	YMC-Pack C8 or Develosil C8-HG-5
December 1, 2023		Assay	YMC-Pack C8 or Develosil C8-HG-5
September 1, 2023	Indocyanine Green	Purity (2) Related substances	YMC-Pack ODS AM
March 1, 2023	Aripiprazole	Purity Related substances, Assay	YMC-Pack Pro C18
Marrala 1, 0000	Febuxostat	Purity Related substances(i), Assay	TSKgel ODS-80Ts
March 1, 2023		Purity Related substances(ii)	Develosil C-30-UG-3
March 1, 2023	Febuxostat Tablets	Identification, Purity Related substances, Uniformity of dosage units, Assay	TSKgel ODS-80Ts
March 1, 2023	Gefitinib Tablets	Assay	Inertsil ODS-3

Pharmaceuticals and Medical Devices Agency Office of Review Management December 1, 2025

Posting Date	Monograph	Test	Brand Name
		Acetic Acid	Spherisorb ODS
March 1, 2023	Goserelin Acetate	Purity	XTerra MS C18
		Assay	XTerra MS C18
March 1 2022	Lawayiaan	Purity Related substances	L-column 2 ODS
March 1, 2023	Lornoxicam	Assay	CAPCELL PAK C18 MGII
March 1, 2023	Lornoxicam Tablets	Purity Related substances, Uniformity of dosage units, Dissolution, Assay	TSKgel ODS-80TM
		Purity(2) Oxalic acid, Purity(3) Related substance A	Hypersil BDS-C18
March 1, 2023	Oxaliplatin	Purity(4) Other related substances, Assay	Hypersil ODS
		Purity(5) Enantiomer	Chiralcel OC-H
March 1, 2023	Oxaliplatin Injection	Purity(1) Related substances, Purity(2) Oxalic acid	Hypersil BDS-C18
		Assay	Hypersil ODS
March 1, 2023		Purity Related substances	TSKgel ODS-100V 3µm
	Tolvaptan	Assay	YMC-Pack ODS-A or YMC-Pack ODS-AM
March 1, 2023	Tolvaptan tablets	Identification, Uniformity of dosage units, Assay	YMC-Pack ODS-A or YMC-Pack ODS-AM

Pharmaceuticals and Medical Devices Agency Office of Review Management December 1, 2025

Posting Date	Monograph	Test	Brand Name
September 1, 2021	Anastrozole Tablets	Uniformity of dosage unit, Dissolution, Assay	Hicrom RPB
September 1, 2021	Bicalutamide Tablets	Assay	Spherisorb ODS2
June 1, 2021	Voglibose Orally Disintegrating Tablets	Uniformity of dosage unit, Dissolution, Assay	YMC-Pack Polyamine II
March 1, 2021	Budesonide	Purity (2) Related substances, Isomer ratio, Assay	Hypersil ODS C18 or Discovery HS C18
		Purity (2) Related substances, Assay	Spherisorb ODS2 5µm 150mm x 4.6mm
September 1, 2020	Temozolomide	Purity (3) Acetonitrile	J&W DB-WAX, 30m x 0.53mm fused silica, 1.0µm film thickness
September 1, 2020	Temozolomide Capsules	Identification, Purity Related substances, Uniformity of dosage units, Assay	Spherisorb ODS2 5µm 150mm x 4.6mm
September 1, 2020	Temozolomide for Injection	Identification, Purity Related substances, Assay	Spherisorb ODS2 5µm 150mm x 4.6mm
June 1, 2020	Anastrozole	Purity (2) Related substances, Assay	Hichrom RPB
March 2, 2020	Oxybutynin Hydrochloride	Purity (2) Related substances	Symmetry C8
September 30, 2019	Cloperastine Fendizoate	Purity (3) 4-Chlorobenzophenone	L-column ODS
September 30, 2019	Cloperastine Fendizoate Tablets	Uniformity of dosage units, Assay	L-column ODS
		Dissolution	TSKgel ODS-80Ts

Pharmaceuticals and Medical Devices Agency Office of Review Management December 1, 2025

Posting Date	Monograph	Test	Brand Name
Santambar 20, 2010	Eribulin Mesilate	Identification (2), Assay (2) Methanesulfonic acid	Zorbax NH2 5μm 4.6mmX25cm
September 30, 2019	Eribulin Meshate	Purity (2) Related substance, Assay (1) Eribulin mesilate	ACE 3 C18 3μm 3.0mmX15cm
September 30, 2019	Phenobarbital Tablets	Uniformity of dosage units, Assay	Mightysil RP-18 GP-Ⅱ
September 30, 2019	Pitavastatin Calcium Orally Disintegrating Tablets	Purity, Uniformity of dosage units, Assay	L-column ODS
	Districegrating rabiets	Dissolution	L-column2 ODS
September 30, 2019	Silodosin Orally Disintegrating Tablets	Identification, Purity Relsted substance, Uniformity of dosage units, Dissolution, Assay	Inertsil ODS-3
September 2, 2019	Bicalutamide	Purity (2) Related substances, Assay	Kromasil C18
September 2, 2019	Celecoxib	Purity (2) Related substances, Assay	SUPELCOSIL LC-DP
September 2, 2019	Dorzolamide Hydrochloride and Timolol Maleate Ophthalmic Solution	Identification (1), Purity (1) Related substances 1, Assay (1) Dorzolamide Hydrochloride	Zorbax Rx−C8 5μm 4.6mm X25cm
		Identification (2), Purity (2) Related substances 2, Assay (2) Timolol Maleate	Inertsil ODS-2 5μm 4.6mmX25cm
September 2, 2019	Ethyl Loflazepate	Purity (4) Related substances	TSKgel ODS-80TM
		Assay	Develosil ODS-7

Pharmaceuticals and Medical Devices Agency Office of Review Management December 1, 2025

Posting Date	Monograph	Test	Brand Name
September 2, 2019	Ethyl Loflazepate Tablets	Unifomity of dosage units, Dissolution, Assay	Inertsil ODS-3
September 2, 2019	Fenofibrate	Purity (2) Related substances, Assay	YMC-Pack Pro C18
		Purity Related substances	CAPCELL PAK C18 UG80
September 2, 2019	Fenofibrate Tablets	Unifomity of dosage units, Dissolution, Assay	YMC-Pack ODS-AM
September 2, 2019	Gefitinib	Purity (2) Related substances, Assay	Inertsil ODS-3 C18
September 2, 2019	Glucagon (Genetical Recombination)	Identification (1)	Lichrospher 100 RP18
September 2, 2019	Glucagon (Genetical Recombination)	Identification (2), Purity, Assay	ACE 3 C18
September 2, 2019	Irinotecan Hydrochloride Injection	Purity Related substances, Assay	Inertsil ODS-2
Santambay 2, 2010	Rilmazafone Hydrochloride Hydrate	Purity (2) Related substances	L-column ODS 5μm 4.6X250 mm
September 2, 2019		Assay	L-column ODS 5μm 4.6X150 mm
September 2, 2019	Rilmazafone Hydrochloride Tablets	Uniformity of dosage units, Dissolution, Assay	L-column ODS 5μm 4.6X150 mm
Santambar 2 2010	Rosuvastatin Calcium	Purity (3) Related substances, Assay	Inertsil ODS-3 C18 3μm 100 Å
September 2, 2019	Rosuvastatiri Galcium	Purity (4) Enantiomer	CHIRALCEL OJ-RH
September 2, 2019	Rosuvastatin Calcium Tablets	Purity Related substances, Assay	Columbus 5μ C18 110A
September 2, 2019		Dissolution	Spherisorb S5 ODS2
L 2 0010	Copovidone	Purity (4)	Inertsil ODS-4
June 3, 2019	Copovidorie	Purity (7)	Inertsil ODS-3
June 3, 2019	Telmisartan and Amlodipine Besylate Tablets	Uniformity of dosage units, Dissolution, Assay	Inertsil C8-3 (5 µ m 3.0 × 75mm)

Pharmaceuticals and Medical Devices Agency Office of Review Management December 1, 2025

Posting Date	Monograph	Test	Brand Name
		Uniformity of dosage units, Assay	YMC-Pack ODS-AM
September 3, 2018	Fludiazepam Tablets	Dissolution	μ BONDASPHERE
		Dissolution	5 μ C18-100A
June 1, 2018	Miglitol Tablets	Uniformity of dosage units, Dissolution,	Asahipak
oune 1, 2016	Migitor rabiets	Assay	NH2P-50 4D
March 1, 2018	Bromfenac Sodium Hydrate	Purity (2) Related substances, Assay	YMC-Pack ODS-A
March 1, 2018	Bromfenac Sodium Ophthalmic Solution	Assay	CAPCELL PAK C18 SG120 Å
March 1, 2018	Cefixime Fine Granules	Purity Related substances, Uniformity	Inertsil ODS-3
March 1, 2018	Cenxime Fine Granules	of dosage unit, Dissolution, Assay	Inertsii ODS-3
	Doripenem Hydrate	Purity (3) Related substances (i),	L-column ODS
Mayala 1 2010		Assay	L-column ODS
March 1, 2018		Purity (3) Related substances (ii)	Kaseisorb LC ODS-SAX Super
		Purity (3) Related substances (iii)	Cadenza CD-C18
	Doripenem for Injection	Purity (2) Related substances (i),	
March 1, 2018		Assay	L-column ODS
		Purity (2) Related substances(ii)	Kaseisorb LC ODS-SAX Super
March 1, 2018	Felodipine	Purity (2) Related substances	L-Column ODS
March 1, 2018	Caladinina Tablata	Uniformity of dosage unit, Dissolution,	L-Column ODS
March 1, 2018	Felodipine Tablets	Assay	
March 1, 2018	Irinotecan Hydrochloride Hydrate	Purity Related substances	Inertsil ODS-2
March 1, 2018	Delegrazina	Purity Related substances, Assay	Mightysil RP18GP
March I, 2018	Polaprezinc		L-Column ODS

Pharmaceuticals and Medical Devices Agency Office of Review Management December 1, 2025

Posting Date	Monograph	Test	Brand Name
March 1, 2018	Polaprezinc Granules	Uniformity of dosage unit, Assay	L-Column ODS Symmetry C18
March 1, 2018	Sodium Valproate Extended-release Tablets A	Uniformity of dosage unit, Dissolution, Assay	Develosil ODS-UG-5
March 1, 2018	Sodium Valproate Extended-release Tablets B	Uniformity of dosage unit, Dissolution, Assay	Develosil ODS-UG-5
March 1, 2018	Telmisartan and Hydrochlorothiazide Tablets	Idendification (1)(2), Uniformity of dosage unit (1)(2), Dissolution (1)(2), Assay (1)(2), Purity Related substances	Inertsil C8-3
March 1, 2018	Valsartan and Hydrochlorothiazide Tablets	Uniformity of dosage unit (1)(2), Dissolution (1)(2), Assay (1)(2)	Nucleosil 100-5 C18
December 1, 2017	Ethylcellulose	Assay	HP-1, Rtx-1
December 1, 2017	Gatifloxacin Hydrate	Purity (3) Related substances, Assay	Inertsil ODS-2
December 1, 2017	Gatifloxacin Ophthalmic Solution	Purity Related substances, Assay	Inertsil ODS-2
December 1, 2017	Hydroxyethylcellulose	Assay	HP-1, Rtx-1
December 1, 2017	Lanoconazole	Purity Related substances, Assay	Inertsil ODS-2
December 1, 2017	Lanoconazole Cream	Assay	Inertsil ODS-2
December 1, 2017	Lanoconazole Cutaneous Solution	Assay	Inertsil ODS-2
December 1, 2017	Lanoconazole Ointment	Assay	Inertsil ODS-2
December 1, 2017	Sitadiatia Phaanhata Hudrata	Purity (1) Related substances, Assay	Supelco Discovery Cyano
December 1, 2017	Sitagliptin Phosphate Hydrate	Purity (2) Enantiomer	Diacel AD-H

Pharmaceuticals and Medical Devices Agency Office of Review Management December 1, 2025

Posting Date	Monograph	Test	Brand Name
December 1, 2017	Sitagliptin Phosphate Tablets	Identification (2), Purity Related substances, Uniformity of dosage unit, Dissolution, Assay	Supelco Discovery Cyano
December 1, 2017	Verapamil Hydrochloride Injection	Assay	Mightysil RP-18GP
September 1, 2017	Cefalotin Sodium for Injection	Purity (2) Related substances, Assay	Inertsil ODS-3
September 1, 2017	Dexamethasone Sodium Phosphate	Purity (5) Related substances, Assay	Waters, Symmetry C18 Column, 100Å
September 1, 2017	Diclofenac Sodium Suppositories	Assay	NUCLEOSIL 5C18 LICHROSPHER PR18-5
September 1, 2017	Ritodrine Hydrochloride Injection	Assay	YMC-Pack Pro C8
September 1, 2016	Acarbose	Purity (3) Related substances, Assay	APS-2 Hypersil
September 1, 2016	Acarbose Tablets	Purity Related substances, Uniformity of dosage units, Assay	APS-2 Hypersil
		Dissolution	Develosil ODS-T
September 1, 2016	Clomipramine Hydrochloride Tablets	Uniformity of dosage units, Assay	Wakosil 10C18
September 1, 2016	Entacapone	Purity (2) Related substances, Assay	XTerra Phenyl, YMC-Pack Ph
September 1, 2016	Entacapone Tablets	Uniformity of dosage units, Assay	XTerra Phenyl
September 1, 2016	Irbesartan Tablets	Uniformity of dosage units, Assay	YMC-Pack ODS-A
September 1, 2016	Pazufloxacin Mesilate Injection	Assay	Nucleosil 100-10C ₁₈
September 1, 2016	Tramadol Hydrochloride	Purity (3) Related substances(ii)	LiChrospher 60 RP-select B

Pharmaceuticals and Medical Devices Agency Office of Review Management December 1, 2025

Posting Date	Monograph	Test	Brand Name
June 1, 2016	Azosemide Tablets	Uniformity of dosage units, Assay	SunFire C18
June 1, 2016	Purified Glucose	Identification (2), Purity (3) Related substances, Assay	BIORAD Aminex HPX-87C
June 1, 2016	Glucose Hydrate	Identification (2), Purity (3) Related substances, Assay	BIORAD Aminex HPX-87C
June 1, 2016	Irbesartan and Amlodipine Besilate Tablets	Dissolution, Identification(1)(2), Uniformity of dosage units(1)(2), Assay(1)(2)	Shim-pack XR-ODS II
March 1, 2016	Pazufloxacin Mesilate	Purity (2) Related substances	Inertsil ODS-2
Iwarch 1, 2010	Fazulioxaciii iviesilate	Assay	Nucleosil 100-10C ₁₈
March 1, 2016	Pentobarbital Calcium Tablets	Uniformity of dosage unit, Assay	L-Column ODS, Inertsil ODS-3
March 1, 2016	Zonisamide	Purity (4) Related substances, Assay	Develosil ODS-5

Disclosure of Information about Columns for Japanese Pharmacopoeia Draft Monographs

Division of Pharmacopoeia and Standards for Drugs,
Office of Standards and Guidelines Development,
Pharmaceuticals and Medical Devices Agency

Please be notified that the division of Pharmacopoeia and Standards for Drugs, Office of Standards and Guidelines Development, Pharmaceuticals and Medical Devices Agency (PMDA) starts to take a measure to in principle disclose the information about columns, such as the name (model number) of columns, for Japanese Pharmacopoeia (JP) draft monographs (hereinafter excluding the draft monographs for crude drugs) under the following rules:

- Concerning tests as a whole using columns in the draft monographs, while further ensuring of transparency of the JP revision process is required, the measure mentioned above is to post the information about the columns that were used for acquisition of the data that were referred in the preparation of the drafts on the PMDA website by PMDA at the start of publication of the relevant drafts for public comments.
- 2. This disclosure is principally provided for ensuring of enhancement of public comments by wide sharing of the information 1 above with stakeholders other than the manufacturers who prepared the drafts at the time of publication for public comments. Thus, addition of information about other alternative columns or update of information associated with technological innovation is in principle not to be done.
- 3. The columns to be disclosed have not been confirmed as applicable to all the samples that could be the subjects to application to the monographs of the relevant drafts.
- 4. The information about columns are released only when cooperation is provided by the manufactures who prepared the drafts.