


Medical Safety Information Pharmaceuticals and Medical Devices Agency http://www.pmda.go.jp/english/service/medical\_info.html

2

(Case 2) The prefilled syringe's gasket was distorted, so the drug solution leaked out from the opening between the syringe barrel and the gasket.

## Precautions when using prefilled syringes (Part 2)

### • When assembling a prefilled syringe, make sure to screw in the plunger by rotating it.




Distortion of gasket may cause drug solution leakage and air bubble entry. When assembling a prefilled syringe, do it correctly following the instructions on the Package Insert, etc.

# Prefilled syringes that are assembled for use.





### Medical Safety Information Pharmaceuticals and Medical Devices Agency http://www.pmda.go.jp/english/service/medical\_info.html






#### Medical Safety Information Pharmaceuticals and Medical Devices Agency http://www.pmda.go.jp/english/service/medical info.html



- We have tried to ensure the accuracy of this information at the time of its compilation but do not guarantee its accuracy in the future.
- This information is not intended to impose constraints on the discretion of healthcare professionals or to impose obligations and responsibility on them, but is provided as a support to promote the safe use of pharmaceuticals and medical devices by healthcare professionals.



Contact: Medical Safety

E-mail: safety.info@pmda.go.jp Information Group http://www.info.pmda.go.jp/